

ABC monoclonal antibody staining of dissected larval and adult CNS

Ito, 20.1.97

For 1st-instar CNS, use PBS instead of PBT throughout the process (small specimen is too fragile for the detergent).

for larva / pupa : Dissect CNS in **PBS**

for adult : Put the animal one minute in **50-70% EtOH** (to wash oil on the body surface)
Transfer to **PBS** and dissect the CNS

Fix: put **100-200** μ l **fixative** into Eppendorf tubes.

most cases: **4% Formaldehyde (FA) / PEM** (100 μ l 37% FA / 900 μ l PEM, 400 μ l 10% FA / 600 μ l PEM)
using "EM grade ultra-pure FA" yields better results

Pool dissected tissues in the **fixative** for up to **10-20 min**

Change to **fresh fixative** and incubate **R.T.** further **50 min** (for 1st instar: 25 min)

Wash: **2 times 5 min** with **PBT** (0.3-0.5 % Triton X /PBS)

For keeping the specimen at 4°C (up to a few days), wash once more with **PBT**.

Block 1: with **0.3% H₂O₂ / Methanol** (100 μ l 30% H₂O₂ / 10 ml Methanol)

R.T. 30 min to kill endogenous peroxidase

Wash: **4 times 5 min** with **PBT**

Block 2: with **1.5% horse serum / PBT** (150 μ l / 10 ml PBT)

R.T. >1 hr to fill nonspecific antibody binding region

1st Ab: **monoclonal antibody** (mouse) (1:500 in 1.5% horse serum / PBT)

4°C Over Night or **R.T. >3 hr**

Retrieve: keep used antibody in 4°C (used antibodies works better; up to a few weeks or a few staining sessions)

Wash: **3 times 5 min** with **PBT**

2nd Ab: **Biotin-conjugated horse anti-mouse IgG** (Vectastain 1:500 in 1.5% horse serum / PBT)

R.T. >2 hr

Retrieve: keep used antibody in 4°C

Wash: **3 times 5 min** with **PBS** (to wash out Triton X before avidin-biotin reaction and DAB staining)

ABC: Mix **10 μ l Solution A**, **10 μ l solution B**, **500 μ l PBS** (Vectastain Elite)

keep **R.T. >30 min** before use

then incubate the specimen **R.T. >1 hr**

Wash: **4 times** with **PBS**

Stain: **0.5-1 mg/ml DAB + 0.03% H₂O₂** (1 μ l 30% H₂O₂ / 1 ml DAB)

1.5% H₂O₂ 1.3 μ l 6.7 μ l

1mg/ml DAB 200 μ l 1000 μ l

R.T. 5-20 min (check the staining under dissecting microscope)

keep hazardous DAB waste !

Wash: **2 times 5 min** with **PBS**

keep hazardous DAB waste !

Post Fix: **2% Glutaraldehyde (GA) / PBS** (80 μ l 25% GA / 920 μ l PBS)

keep hazardous DAB waste !

R.T. > 1 hr

Wash: **3 times 5 min** with **H₂O**

keep hazardous DAB waste !

Dehydrate: **30%, 50%, 70%, 90%** each **3 min** -> **100% 5 min** -> **100% 10 min** -> **Xylen 15 min** -> **Araldit 30 min**

change **Araldit** again

Mount: in **Araldit**

Sigma DAB Tetrahydrochloride Tablet (D5905, 100 tablets)

antibody / Fluorescent staining of dissected larval and adult CNS

Ito, 21.1.99

For 1st-instar CNS, use PBS instead of PBT throughout the process (small specimen is too fragile for detergent).

for Larva : Dissect CNS in **PBS**

for adult : Put the animal one minute in **50-70% EtOH** (to wash oil on the body surface)

Transfer to **PBS** and dissect the CNS

Fix: **4% Formaldehyde (FA) / PEM** (100µl 37% FA / 900µl PEM, 400 µl 10% FA / 600µl PEM)
R.T. 50 min (1st instar: 25 min) *with gentle agitation*

Wash: **4 times 5 min** with **PBT**

Block: with **10% goat serum / PBT**

R.T. >1 hr to fill nonspecific antibody binding region

1st Ab: **antibody in 10% goat serum / PBT**

4°C Over Night or **R.T. >3 hr**

Retrieve: keep used antibody in 4°C (used antibodies works better; up to a few weeks or a few staining sessions)

Wash: **3 times 5 min** with **PBT**

2nd Ab: **Cy3-conjugated anti IgG** (Jackson, **1:100** in **10% goat serum / PBT**)

(anti mouse IgG for monoclonal, anti rabbit IgG for most polyclonal)

R.T. >2 hr

during incubation, put cover onto the Eppendorf tubes to keep specimen in the dark (in all the steps)

Wash: **3 times 5 min** with **PBS**

Dehydrate: in **50% Glycerol / PBS** for 1 hour

in **50% antifade Glycerol / PBS** for 30 min

Mount: in **50% antifade Glycerol / PBS**

Store: keep in refrigerator until observation (up to 1 week)

antifade 0.25 % n-propyl gallate 50% Glycerol, PBS
dissolve 25 mg n-propyl gallate in 5 ml PBS
add 5 ml Glycerol
adjust pH to 8.0-8.6 with 1 N NaOH (50-80 µl)

1st : Cappel Rabbit polyclonal anti β-gal
Promega Mouse monoclonal anti β-gal
Sigma Mouse monoclonal anti bovine-tau
Sigma Rabbit polyclonal anti-WGA

Cosmo Bio: 55693 (3712-0081)
Z 3781
T5530
T4144

2nd: Jackson AffinitiPure Goat Anti-Mouse IgG (H+L) Cy3-conjugated 115-165-146
Jackson AffinitiPure Goat Anti-Rabbit IgG (H+L) Cy3-conjugated 115-165-144

(1612-0081) Cappel: affinity purified goat anti rabbit IgG (H+L) - FITC conjugated (whole) Cosmo Bio: 55662

Zymed: affinity purified goat anti rabbit IgG (H+L) - FITC conjugated (whole) Cosmo Bio: 62-6111

Preabsorbition: If background is high, incubate the diluted antibody for >2 hrs with fixed embryos (dechorionate with breach for 4 min, fix with heptan / 10% FA for 30 min , deviterinize with heptan/methanol, wash with methanol and then with PBS)

GFP preparation of dissected larval and adult CNS

Ito, 20.1.99

- for Larva : Dissect CNS in **PBS**
for adult : Put the animal one minute in **50-70% EtOH** (to wash oil on the body surface)
Transfer to **PBS** and dissect the CNS
- Fix: put **100-200 µl 4% Formaldehyde (FA) / PEM** into Eppendorf tubes.
(100µl 37% FA / 900µl PEM, or 400 µl 10% FA / 600µl PEM)
Pool dissected tissues in the **fixative** for up to **10-20 min**
Change the **fresh fixative** and incubate **R.T.** further for **50 min** *with gentle agitation*
during incubation, put cover onto the Eppendorf tubes to keep specimen in the dark (in all the steps)
- Wash: **3 times 5 min** with **PBS**
- Dehydrate: in **50% Glycerol / PBS** for 1 hour
in **50% antifade Glycerol / PBS** for 30 min
- Mount: in **50% antifade Glycerol / PBS**
- Store: keep in refrigerator until observation (up to 1 week)
- antifade** 0.25 % n-propyl gallate 50% Glycerol, PBS
dissolve 25 mg n-propyl gallate in 5 ml PBS
add 5 ml Glycerol
adjust pH to 8.0-8.6 with 1 N NaOH (50-80 µl)

Anti β -gal DAB staining of dissected larval and adult CNS

Ito, 4.6.95

For 1st-instar CNS, use PBS instead of PBT throughout the process (small specimen is too fragile for detergent).

- for Larva : Dissect CNS in **PBS**
for adult : Put the animal one minute in **50-70% EtOH** (to wash oil on the body surface)
Transfer to **PBS** and dissect the CNS
- Fix: put **100-200 μ l fixative** into Eppendorf tubes.
most cases: **4% Formaldehyde (FA) / PEM** (100 μ l 37% FA / 900 μ l PEM, 400 μ l 10% FA / 600 μ l PEM)
some cases: **1% Glutaraldehyde (GA) / PBS** (40 μ l 25% GA / 960 μ l PBS)
GA works well with anti β -gal, but some antibodies don't work with GA at all.
Pool dissected tissues in the **fixative** up to **10-20 min**
Change the **fresh fixative** and incubate **R.T.** further
FA: **50 min**, FA for 1st instar: 25 min, GA: 10 min
- Wash: **2 times 5 min** with **PBT** (0.3-0.5 % Triton X /PBS)
For keeping the specimen at 4°C (up to a few days), wash once more with **PBT**.
- Block 1: with **0.3% H₂O₂ / Methanol** (1/100 volume of 30% H₂O₂)
R.T. 30 min to kill endogenous peroxidase
- Wash: **4 times 5 min** with **PBT**
- Block 2: with **10% goat serum / PBT**
R.T. >1 hr to fill nonspecific antibody binding region
- 1st Ab: **rabbit polyclonal anti β -gal** (Cappel 1:1000 in 10% goat serum / PBT)
4°C Over Night or **R.T. >4 hr**
- Retrieve: keep used antibody in 4°C (used antibodies works better; up to a few weeks or a few staining sessions)
Wash: **3 times 5 min** with **PBT**
- 2nd Ab: for HRP: **Peroxidase-conjugated anti rabbit IgG** (Cappel, 1:500 in 10% goat serum / PBT)
for ABC: **Biotin-conjugated anti rabbit IgG** (Vectastain 1:500 in 10% goat serum / PBT)
R.T. >2 hr
- Retrieve: keep used antibody in 4°C
Wash: **3 times 5 min** with **PBS** (to wash out Triton X before avidin-biotin reaction and DAB staining)
- HRP: goto **Stain**
ABC: Mix **10 μ l Solution A, 10 μ l solution B, 500 μ l PBS** (Vectastain Elite)
keep **R.T. >30 min** before use
then incubate the specimen **R.T. >1 hr**
Wash: **4 times** with **PBS**
- Stain: **0.5-1 mg/ml DAB + 0.03% H₂O₂** (1 μ l 30% H₂O₂ / 1 ml DAB)
1.5% H₂O₂ 1.3 μ l 6.7 μ l
1mg/ml DAB 200 μ l 1000 μ l
R.T. 5-20 min (check the staining under dissecting microscope) **keep hazardous DAB waste !**
- Wash: once **5 min** with **PBS** **keep hazardous DAB waste !**
- Post Fix: **2% Glutaraldehyde (GA) / PBS** (80 μ l 25% GA / 920 μ l PBS) **keep hazardous DAB waste !**
R.T. > 1 hr
- Wash: **3 times 5 min** with **H₂O** **keep hazardous DAB waste !**
- Dehydrate: in **50% Glycerol / H₂O** for 1 hour
in **80% Glycerol / H₂O** for 30 min
- Mount: in **80% Glycerol / H₂O**
- or:
For plastic mount:
Dehydrate: **30%, 50%, 70%, 90%** each **3 min** -> **100% 5 min** -> **100% 10 min** -> **Xylen 15 min** ->
Mount: in **Araldit**
- Stock
1st : Cappel Rabbit polyclonal anti β -gal Cosmo Bio: 55693 (3712-0081)
2nd affinity purified goat anti rabbit IgG (H+L) (whole)
-- peroxidase conjugated : Cappel : Cosmo Bio: 55689 (3612-0081), Zymed : Cosmo Bio: 62-6120
-- biotinilated (whole) Cappel : Cosmo Bio: 55699 (8612-3711), Vector : Vectastain Elite
- DAB tablet: Sigma DAB Tetrahydrochloride Tablet Sigma D-5905
Goat Serum: CMN Chemicon Cosmo Bio S-26 017

Antibody staining of Embryos (with Vectastain kit)

ITO Kei, 1995 April

- Dechorionate: collect eggs with a nylon mesh attached to the bottom of a funnel
soak the eggs in **7.5% Bleach in H₂O** (50% Chrolox, 50% H₂O)
R.T. 4 min
- Wash: with **H₂O**
- Fix: Transfer into an Eppendorf tube containing
FA / Heptan (1 ml)
400 µl PBS, 100µl 37 % Formaldehyde, 500 µl n-Heptan
(or 100 µl PBS, 400 µl 10 % EM grade Formaldehyde, 500 µl n-Heptan)
R.T. 30-40 min
- Deviterinize: remove the FA+PBS phase (lower part)
add the same amount of **100% Methanol**
shake well by hand or by Vortex
remove the upper phase (Methanol)
add another volume of **Methanol**
shake well by hand or by Vortex
until the majority of the eggs go down to the bottom
remove heptan and Methanol solutions
- Wash: **2 times with Methanol**
until eggs become non-sticky
Keep embryos at this stage at 4°C (up to a few weeks)
- Before staining, wash **4 times with PBT** (0.3-0.5% Triton X-100 in PBS)
- Block: with **10% goat serum / PBT**
R.T. >1 hr to fill nonspecific antibody binding region
- 1st Ab: **monoclonal antibody** (mouse) or **polyclonal antibody** (rat, rabbit, goat etc.)
(in **10% goat serum / PBT**) *preabsorbed!*
4°C Over Night or **R.T. >3 hr**
- Retrieve: keep used antibody in 4°C (used antibodies works better; up to a few weeks or a few staining sessions)
- Wash: **3 times with PBT**
- 2nd Ab: Biotin-conjugated **antibody (1:500 in 10% goat serum / PBT; Vectastain)**
for monoclonal 1st: anti mouse Ig, for polyclonal 1st: anti rat, rabbit, goat,, Ig
preabsorbed!
R.T. >2 hr
- Retrieve: keep used antibody in 4°C
- Wash: **3 times with PBS** (to wash out Triton X before avidin-biotin reaction)
- ABC: Mix **10µl Solution A, 10µl solution B, 500µl PBS** (Vectastain Elite)
keep **R.T. >30 min** before use
then incubate the specimen **R.T. >1 hr**
- Wash: **4 times with PBS**
- Stain: **0.5-1mg/ml DAB + 0.01% H₂O₂**
1.5% H₂O₂ 1.3µl 6.7µl
1mg/ml DAB 200µl 1000µl
R.T. 5-20 min (check the staining under dissecting microscope) ***keep hazardous DAB waste !***
- Wash: **2 times with PBS** ***keep hazardous DAB waste !***
- Wash: once with **H₂O** ***keep hazardous DAB waste !***
- Dehydrate: in **50% Glycerol / H₂O** for 1 hour
in **80% Glycerol / H₂O** for 30 min
- Mount: in **80% Glycerol / H₂O**
- For plastic mount:
- Dehydrate: **30%, 50%, 70%, 90%** each **3 min** -> **100% 5 min** -> **100% 10 min** -> **Xylen 15 min** ->
- Mount: in **Araldit**

Preabsorbtion: incubate the diluted antibody for >2 hrs with fixed embryos (dechorionate with breach for 4 min, fix with heptan / 10% FA for 30 min , deviterinize with heptan/methanol, wash with methanol and then with PBS)

X-gal staining of embryo

Ito, 4.6.95

- Dechorionate: put the eggs in **7.5% Bleach** (in **H₂O**)
R.T. 4 min
- Wash: collect the eggs with a nylon mesh
with **H₂O**
transfer to an Eppendorf tube
- Fix: **GDA / Heptan** (300-500 µl)
45 ml PBS, 5 ml 25% Glutardiardehid, 50 ml n-Heptan
store at 4°C, shake well just before use and take the upper phase
R.T. 15 min
- Deviterinize: (divitelinize by hand works better)
remove **GDA / PBS** phase (lower part)
add the same amount of **80% EtOH** (in H₂O)
shake well by hand or by Vortex
(Although Ethanol deviterinization is not so effective as with
Methanol, the loss of β-galactosidase activity is smaller.)
remove the upper phase (EtOH)
add another **80% EtOH**
shake well by hand
until the majority of the eggs go down to the bottom
remove heptan and EtOH solutions
- Wash: several times with **80% EtOH**
until eggs become non-sticky
wash 2 times with **PBT** (0.5% Triton X-100 in PBS)
- Stain: heat **X-gal buffer** to **65°C**
buffer: 100 ml 0.1 M Na phosphate buffer
860 mg NaCl
20.3 mg MgCl₂
109 mg K₃[Fe(III)CN₆]
140 mg K₄[Fe(II)CN₆]
store at 4°C
add **1/100** volume of **20% X-gal**
200 mg X-gal in 1 ml DMF, store at 4°C
mix well and *gradually* cool down to room temperature
If precipitation appears, heat and mix the solution again.
apply the staining solution to the eggs
37°C 2 hr to Over Night (with gentle shaking)
- Wash: **3 times 5 min** with **PBS**
- Dehydrate: in **50% Glycerol / H₂O** for 1 hour
in **80% Glycerol / H₂O** for 30 min
- Mount: in **80% Glycerol / H₂O**
- or:
For plastic mount:
EtOH Dehydrate: **70%, 90%** each 5 min -> **100%** 5 min -> **100%** 10 min -> **Xylen** 15 min
Mount: in **Araldit**

X-gal staining of dissected larval and adult CNS

Ito, 4.6.95

- for Larva : Dissect CNS in **PBS**
for adult : Put the animal **one** minute in **50-70% EtOH** (to wash oil on the body surface)
Transfer to **PBS** and dissect the CNS
- Fix: put **100-200** μ l **fixative** into Eppendorf tubes.
1% Glutaraldehyde (GA) / PBS (40 μ l 25% GA / 960 μ l PBS)
Pool dissected tissues in the **fixative** up to **10-20 min**
Change the **fresh fixative** and incubate **R.T. 15 min**
- Wash: **2 times 5 min** with **PBT** (0.3-0.5 % Triton X /PBS)
- Stain: heat **X-gal buffer to 65°C**
buffer: 100 ml 0.1 M Na phosphate buffer
860 mg NaCl
20.3 mg MgCl₂
109 mg K₃[Fe(III)CN₆]
140 mg K₄[Fe(II)CN₆]
- store at 4°C
add **1/400** volume of **20% X-gal**
200 mg X-gal in 1 ml DMF, store at 4°C
mix well and *gradually* cool down to room temperature
If precipitation appears, heat and mix the solution again.
- apply the staining solution to the eggs
37°C 2 hr to Over Night (with gentle shaking)
- Wash: **3 times 5 min** with **PBS**
- Dehydrate: in **50% Glycerol / H₂O** for 1 hour
in **80% Glycerol / H₂O** for 30 min
- Mount: in **80% Glycerol / H₂O**
- or:
For plastic mount:
EtOH Dehydrate: **70%, 90%** each 5 min -> **100%** 5 min -> **100%** 10 min -> **Xylen** 15 min
Mount: in **Araldit**

Fixative

4 % Formaldehyde (FA) good for the first trial

10% Formalin / PBS
 PBS Tablets pH7.4 (25 degree) (Sigma P4417)

4 % Ultra pure Formaldehyde (FA) usually the good and most convenient

40% EM grade Ultra pure Formaldehyde / PEM
 Formaldehyde, Methanol Free, 10%, Ultra pure EM grade (Polyscience, 04-1040-18)

4 % Paraformaldehyde (PFA) good for fixing protein

mix: paraformaldehyde	4		g 1		g
0.1 M Phosphate Buffer	90		ml 22		ml
warm to 60°C (ca. 20 min.)					
cool and mess up to	100		ml 25		ml
store at 4°C (<i>max 1 week</i>)					

PLP (Periodate–Lysine–Paraformaldehyde)

PFA for protein, periodate-lysine for sugar

Stock A:

mix: Lysine–HCl	1.827 g		0.365 g		
H ₂ O 50	ml		10 ml		
add 0.1 M Na ₂ HPO ₄ to pH 7.4					
mess up with 0.1 M Phosphate Buffer to	100		ml 20		ml
store at 4°C (<i>max 10 days</i>)					

Stock B:

mix: Paraformaldehyde (PFA)	8		g 0.2 g		
H ₂ O 100	ml		5 ml		
warm to 60°C (ca. 20 min.)					
add 1N NaOH			a few drops to dissolve PFA		
filtrate					
store at 4°C (<i>max only a few days</i>)					

Just before use:

mix Stock A	7.5 ml	9 ml	3 ml
Stock B	2.5 ml	3 ml	1 ml
NaIO ₄	21.4 mg	25.68 mg	8.56 mg
total	10 ml	12 ml	4 ml
<i>(final 0.01 M NaIO₄ – 0.075 M Phosphate – 2 % PFA pH. 6.2)</i>			

0.1 M Na-Phosphate Buffer


	for 1x soln.	for 10x soln.
mix Na ₂ HPO ₄ · H ₂ O	10.22 g	51.1 g
or Na ₂ HPO ₄ · 2H ₂ O	12.82 g	64.1 g
or Na ₂ HPO ₄ · 7H ₂ O	19.30 g	96.5 g
or Na ₂ HPO ₄ · 10H ₂ O	25.79 g	128.95 g
AND NaH ₂ PO ₄ · H ₂ O	3.86 g	19.3 g
or NaH ₂ PO ₄ · 2H ₂ O	4.37 g	21.85 g
total	1 l	500 ml

PEM Buffer

	for 100 ml	for 500 ml
mix 0.1 M PIPES (302.37)	3.0237 g	15.1185 g
2 mM EGTA (380.35)	76.07 mg	380.35 mg
1 mM MgSO ₄ · 7H ₂ O (246.48)	24.648 mg	123.24 mg
adjust pH to 6.95 with NaOH (or HCl)		

Preparation of Araldit

Araldit CY 212	(Serva, #13824, 1 kg)
Araldit hardener HY 964	(Serva, #13826, 1 kg)
Araldit softener (Phtalic Acid Dibuthyl Ester: PADE)	(Serva, #32805, 250 ml)
Araldit Beschleuniger DY 964	(Serva, #13825, 100 g)



For semi-thin sectioning

100 ml Araldit CY 212 into a 500 ml glass cup.
Add **100 ml Araldit hardener HY 964** slowly.
Put the rotor of a low-speed electric mixer in it.
Let the mixer slowly rotate for at least **2 hrs.** (avoid making bubbles in the solution.)
Add **10 ml Softener PADE** very slowly (drop by drop) with the mixer rotating.
Adding as little as 500 ul at once will result in uneven solution,
which is not usable any more.
Mix for another **15 min.**
Add **5 ml Beschleuniger DY 964** very slowly (drop by drop).
Mix for another **30 min.**

Aliquot the Araldit in 10-20 ml disposable syringe.
Store frozen at -20°C.

For whole-mount preparation

Softener PADE is not needed. Mix CY 212, HY 964, and DY 964 in this order.
Araldit for semi-thin sectioning is also applicable to whole-mount preparation, although it costs a bit more.

Araldit monomer is harmful. Use gloves while handling. Wipe out with acetone in case of an accident.
Beschleuniger is a highly volatile material; work in a draft chamber.

After use, wipe the glass cup and mixing rotor with acetone and keep them at 63°C over night until remaining araldit become polymerized.

Making Semi-thin Sections

Ito, 4.6.95

Larva, Pupa, Imago:

Dissect animals in Ringer or PBS

Embryo:

Prefix with 25% GDA / Heptan

R.T. 10 min

Wash 3 times with PBS (Remove Heptan)

Fixation 1:

Fix with 6% GDA (Glutardialdehyd)/PBS

4°C, 1-4 hr

Wash 2 times with PBS

Fixation 2:

Fix with 2% OsO₄/PBS

4°C, in the Dark, Over Night

CAUTION: AVOID SKIN and GAS CONTACT

Wash 3 times with PBS

Dehydrate and Embed:

30%, 50%, 70%, 90% each 3 min -> 100% 5 min -> 100% 10 min ->

Xylen 10 min 2 times ->

Araldit/Xylen (1:1) 6hr-Over Night in the Abzug (let Xylen evaporate) ->

Araldit 8-12 hr (let residual Xylen evaporate) ->

Araldit in the mold

Polymerize 45°C, 12 hr -> 65°C, 36-48 hr

Section

Staining

Richardson's Solution

1% Toluidine blue and 1% azur II - 1% Borax

60°C, 15-60 sec.

or

TB Solution

1 ml 0.2% Toluidine blue / water

1 ml 1% Methylen blue / 1% Borax

18 ml water

60°C, 30 sec. - 2 min

Wash with running tap water

Wash with distilled water

Dry

Mount with Araldit

Anti β -gal staining of larval CNS original in Mainz

Dissect CNS

Fix: 1% Glutardialdehyde (GDA) / PBS (40 μ l 25% GDA / 960 μ l PBS)
R.T. 10 min
GDA does not work with other antibodies! Anti β -gal works with Formaldehyde.
Use Formaldehyde (FA) when staining embryos.

Wash: 2 times with PBS

Block 1: with 0.3% H₂O₂ / Methanol
R.T. 30 min to kill endogenous peroxidase

Wash: 3 times with PBS

Block 2: with 10% lamb serum / PBS
R.T. >1 hr to fill nonspecific antibody binding region

1st Ab: anti β -gal (rabbit polyclonal) (1:7000 in 10% lamb serum / PBS)
1:14 preabsorbed antibody solution (1:500)
4°C Over Night

Retrieve: keep used antibody in 4°C

Wash: 3 times with PBS

2nd Ab: Biotin conjugated anti rabbit IgG (1:500 in 10% lamb serum / PBS)
(Vector stain) preabsorbed!
R.T. 2 hr

Retrieve: keep used antibody in 4°C

Wash: 3 times with PBS

ABC: Mix 10 μ l Solution A, 10 μ l solution B, 500 μ l PBS
keep R.T. >30 min before use
then R.T. 1 hr with the specimen

Wash: 4 times with PBS

Stain: 1mg/ml DAB + 0.01% H₂O₂
1.5% H₂O₂ 1.3 μ l 6.7 μ l
1mg/ml DAB 200 μ l 1000 μ l
R.T. 5-10 min

Wash: 2 times with PBS

Post Fix: 1% Glutardialdehyde (GDA) / PBS (40 μ l 25% GDA / 960 μ l PBS)
R.T. 30 min

Wash: once with Water

Dehydrate: 30%, 50%, 70%, 90% each 3 min -> 100% 5 min
-> 100% 10 min -> Xylen 15 min -> 50% Araldit / Xylen O/N
-> Araldit 6 hr -> Araldit